



Hinnant Prosthetics

Experience Our Road to Prosthetic Excellence

120 East Kingston Avenue, Charlotte, NC 28203 Ph: 407/375/2587 Fax: 704/333/4429

Sock Ply Combination Chart

1 ply	= 1 ply
1 ply + 1 ply	= 2 ply
3 ply sock	= 3 ply
3 ply + 1 ply	= 4 ply
3 ply + 1 ply + 1 ply	= 5 ply
5 ply sock	= 5 ply
5 ply + 1 ply or (3ply + 3 ply)	= 6 ply
5 ply + 1 ply + 1 ply or (3ply + 3 ply + 1 ply)	= 7 ply
5 ply + 3 ply	= 8 ply
5 ply + 3 ply + 1 ply or (3 ply + 3 ply + 3 ply)	= 9 ply
5ply + 5 ply or (3ply + 3 ply + 1 ply)	= 10 ply
5 ply + 5 ply + 1 ply	= 11 ply
5 ply + 5 ply + 1 ply + 1 ply	= 12 ply
5 ply + 5 ply + 3 ply	= 13 ply
5 ply + 5 ply + 3 ply + 1 ply	= 14 ply
5 ply + 5 ply + 5 ply	= 15 ply
5 ply + 5 ply + 5 ply + 1 ply	= 16 ply
5 ply + 5 ply + 5 ply + 1 ply + 1 ply	= 17 ply
5 ply + 5 ply + 5 ply + 3 ply	= 18 ply
5 ply + 5 ply + 5 ply + 3 ply + 1 ply	= 19 ply
5 ply + 5 ply + 5 ply + 5 ply	= 20 ply

ETC.

Sock Band Colors:

1 ply = White
3 ply = Yellow
5 ply = Green

Should your socks have a number on them and not a colored band:

1 = 1 ply
3 = 3 ply
5 = 5 ply

<http://www.hinnantprosthetics.com>